



# **An overview of isothermal amplification platforms in general and AmplifyRP® specifically**

*Paul Russell, Nathan McOwen, Shulu Zhang,  
Bryant Davenport, and Rugang Li*

# Why Isothermal?

- Molecular detection is all the rage
  - “I’ve got a PCR for that!”
  - Lateral flow or ELISA are not always enough
  - Confirmation of immunological results
  - Sometimes only game in town
    - Hybridization, PCR and qPCR are lab based methods
      - Labor intensive, time consuming, and expensive
  - Easy, fast, cheap (relatively)
  - Potential for being a Point of Care technology

# Various isothermal detection platforms

**Table 1** Various isothermal detection platforms and components of the assay that are amplified

Section	Platform	Amplified component	Amplification catalyst
3.01	NASBA	Complementary sequence of target (RNA)	Enzymatic
3.02	SDA	Probe	Enzymatic
3.03	LAMP	Probe	Enzymatic
3.04	Invader assay	Signal	Enzymatic
3.05	RCA	Probe	Enzymatic
3.06	SMART	New RNA	Enzymatic
3.07	HDA	Target and the complementary sequence	Enzymatic
3.08	RPA	Target and the complementary sequence	Enzymatic
3.09	NESA, NEANA	Signal	Enzymatic
3.10	Exo III aided target recycling	Signal	Enzymatic
3.11	Junction or Y probes	Signal	Enzymatic
3.12	Reactivation of enzymatic activity	Signal	Enzymatic
3.13	Ultrasensitive miRNA detection	Probe and signal	Enzymatic
3.14	Split DNAzyme	Signal	Non-enzymatic
3.15	RNA cleaving deoxyribozymes (CMB, TASC)	Signal	Non-enzymatic
3.16	Template-directed chemical reaction	Signal	Non-enzymatic
3.17	Non-covalent DNA catalytic reaction	Signal	Non-enzymatic
3.18	HCR	Signal	Non-enzymatic
3.19	Self-assembled DNA	Signal	Non-enzymatic

# RPA Method of Detection - AmplifyRP™


- Based upon the recombinase-polymerase methodology
  - Amplify DNA at a constant temperature (e.g. 37-39°C)
- Uses a mix of enzymes/reagents in a lyophilized pellet
  - Recombinases, DNA polymerase, SSBP and other proteins
- Rapid detection of DNA or RNA
  - Positive results can be seen between 5-20 minutes, depending on detection method

# AmplifyRP Basics

## The RPA Cycle



All steps operate at low constant temperature (optimum 37°C)

a. Recombinase / oligonucleotide primer complexes form and target homologous DNA



Citation: Piepenburg O, Williams CH, Stemple DL, Armes NA (2006) DNA detection using recombination proteins. PLoS Biol 4(7): e204. DOI: 10.1371/journal.pbio.004020

## Not Just the Basics – Two Formats – They Differ How?



**Acceler8™**  
**Think Nested PCR**

● Fluorophore    ^ Abasic site    ■ 3' block    nfo




# What does it look like?



## Procedure



- Extract sample (1:10 – 1:20 w/v) using sample extraction buffer
- Mesh bags, tube and pestle, bead beater, or any method will work





- Rehydrate reaction pellet






- Transfer 1ul of crude sample extract to the rehydrated reaction pellet
- Mix



# Detection of AmplifyRP XRT - Instrument Agnostic



# Acceler8™ is endpoint detection by ImmunoStrip®




**Detection Result**

- How it works:
  - Uses a labelled streptavidin to capture the biotinylated 3' primer and monoclonal anti-dye to capture FITC 5' primer of the amplimer.


# AmplifyRP Hybrid Assays

Same result on DNA dilution series: 100fg – 1pg sensitivity

Real-time fluorescent detection




Endpoint Immunostrip detection



# Sensitivity: Grapevine Red Blotch AmplifyRP vs. PCR

AmplifyRP Acceler8 is at least 100X more sensitive than PCR when testing purified GRBaV-infected grapevine leaf DNA





# ***Candidatus liberibacter asiaticus* (Las)**

- Citrus Greening (HLB)
  - Kills citrus trees.
  - Asia, Africa, S. America (Brazil), Mexico, USA (Florida), Europe (?).
  - 3 strains Las, Lam, Laf.
  - Insect vector (psyllid).



# Assay Development

- Product definition
  - Specificity, sensitivity, cross-reactivity, test format.
  - Other tests on the market.
  - Collaborators and sources of materials (development and validation).
- Primer/probe design
  - Sequence data, lots of it.
  - Lineups, blasts, choice of target region.
  - Choose (no algorithm) and test combinations,
  - Optimize.
- Produce pilot lots and validate.



(repeat).


# Las Acceler8™ Cross-reactivity



Las cross-reactivity with HLB species and *Xanthomonas citri* strains.

Left to right: Las DNA positive control, host DNA negative control, Lam sample, Laf sample, Lsol sample, *X. citri* strain 1, *X. citri* strain 2, *X. citri* strain 3.

# Comparison of Acceler8™ and qPCR



\* Leaf and psyllid samples, DNA, and qPCR data courtesy of Dr. Mike Irey, US Sugar

# Las Acceler8™ Real World Samples

Template		HLB Status			
		Positive (Cq range)		Negative (Cq range)	
		qPCR	Acceler8™	qPCR	Acceler8™
<b>Trial A</b>	<b>DNA</b>	18 (22.49 – 27.16)	18 (NA)	30 (40)	30 (NA)
<b>Trial B</b>	<b>DNA</b>	76 (21.52 – 30.1)	76 (NA)	18 (40)	18 (NA)
	<b>Crude Extract</b>	(NA)	76	(NA)	18
<b>Trial C</b>	<b>DNA</b>	16 (25.22 – 33.0)	16 (NA)	78 (40)	78 (NA)
<b>Trial D</b>	<b>Crude Extract</b>	ND*	31%	ND*	69%

\* qPCR not done. Previous year's data gave an expectation of 20% positives, but is not confirmed this year.

# Summary

- The test is robust: Works very well with extremely variable samples and targets using a crude extract
- Fast: 30 - 45 min start to finish, including sample prep
- Can be done on-site with ImmunoStrip® detection for a visual read or a portable reader for real-time results
- Specific for *L. asiaticus*; no cross-reactivity with other citrus pathogens
- Very sensitive assay: equivalent to qPCR with real world plant tissue and psyllid samples
  - Sensitivity detecting Purified PCR fragment = Approximately 48 copies



# Research and Products

## Pathogen

## Format

\* RNA test

*Phytophthora ramorum*

[XRT] [Acceler8™]

*Phytophthora kernoviae*

[XRT] [Acceler8™]

*Phytophthora palmivora*

[XRT] [Acceler8™]

*Fusarium*, Race 4

[Acceler8™]

*Candidatus liberibacter asaiicus*

[XRT] [Acceler8™]

*Candidatus liberibacter solanacearum*

[Acceler8™]

Phytoplasma [aster yellows]

[XRT] [Acceler8™]

*Clavibacter michiganensis* subsp. *michiganensis*

[XRT]

*Clavibacter michiganensis* subsp. *sepidonicus*

[XRT] [Acceler8™]

*Pseudomonas syringae* pv. *tomato*

[XRT] [Acceler8™]

*Xanthomonas campestris* pv. *vesicatoria* (Xcv)

[XRT] [Acceler8™]

Banana bunchy top virus

[XRT] [Acceler8™]

Plumbox virus\*

[XRT] [Acceler8™]

Little cherry virus 2\*

[Acceler8™]

Tomato chlorotic dwarf viroid\*

[XRT] [Acceler8™]

Potato spindle tuber viroid\*

[XRT] [Acceler8™]

Discovery Format

[XRT] [Acceler8™]